Pre-Mixed pH Calibration

Solutions

PHA Series Basic Unit

- ✓ Pre-Mixed pH Solutions ✓ Two Bottle Sizes **Available**
- Formulated to provide precise pH determinations, these buffer solutions are ideal for calibrating electrodes and pH meters. The buffers are colour-coded to avoid error and accurate to within±0.02 pH at 25°C. Each buffer contains a preservative/ mould inhibitor and a shrink-sealed cap to prevent contamination and leakage. Temperature compensation charts are included with each buffer solution.



IN STOCK FOR FAST DELIVERY!

To Order (Specify Model Number)			
Model No.	Price	Description	
PHA-4	£4.10	4.01 pH buffer solution 475 ml bottle	
PHA-4-GAL	14.75	4.01 pH buffer solution 3.8 L bottle	
PHA-7	4.10	7.00 pH buffer solution 475 ml bottle	
PHA-7-GAL	14.75	7.00 pH buffer solution 3.8 L bottle	
PHA-10	4.10	10.01 pH buffer solution 475 ml bottle	
PHA-10-GAL	14.75	10.01 pH buffer solution 3.8 L bottle	

For Traceable certificate of conformance, contact sales and add £20.50 to price Ordering Example: PHA-7-GAL, 7.00 pH buffer solution 3.8 L bottle, £14.75

ating Buffer Calibrating Buffer Ready to use. ordent: Potasium Hydrogen Plateina, CASe 24-7 in Desembed Water CASe 7732-18-5 se ordentroped added. Just Consones - 20 mt.

PHA-20 Series

✓ Single-Use Packages

The PHA-20 pH buffer sachets are ideal for field use in calibration of pH electrodes and instrumentation. The sealed foil package prevents light, air and other contaminants from affecting the buffer solution. Each package contains enough solution for one calibration. Each box contains 20 sachets.

IN STOCK FOR FAST DELIVERY!

To Order (Specify Model No.)				
Model No.	Price	Description		
PHA-20-4	£14.75	pH buffer 4, 20 sealed pouches		
PHA-20-7	14.75	pH buffer 7, 20 sealed pouches		
PHA-20-10	14.75	pH buffer 10, 20 sealed pouches		

Ordering Example: PHA-20-10, pH buffer solution for pH 10, 20 sealed sachets. £14.75